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Patient Safety Resources

Patient Safety Resources

Novel Coronavirus (COVID-19) Resource Center

Perioperative Considerations for the 2019 Novel Coronavirus (COVID-19)

FAQ on Anesthesia Machine Use, Protection, and Decontamination During the COVID-19 Pandemic

RAPID Response to questions from readers (formerly Dear SIRS)

Current Anesthesia Drug Shortages

Perioperative Multi-Center Handoff Collaborative

Clinical Safety Tools

FAQ ON ANESTHESIA MACHINE USE, PROTECTION, AND DECONTAMINATION DURING THE COVID-19 PANDEMIC

Last updated: March 23, 2020

Contributors: Jeffrey Feldman, MD, Professor of Anesthesiology and APSF Chair of Committee on Technology; Robert "Butch" Loeb, MD, Professor of Anesthesiology University FL; James Phillip, MD, Professor of Anesthesiology Brigham & Women's

Disclaimer: We aim to present our recommendations and a forum for questions. Given the novelty of COVID-19, best-available clinical evidence is supported from anecdotal reports from China, South Korea, Italy and previous epidemics like SARS and MERS. The APSF does not support or endorse any specific product, equipment or trademarked technique. We strongly promote consistency with your governing bodies and organizations such as the CDC, WHO, ASA, AANA, and AAAA.

The following FAQ information is intended to provide the best guidance possible for safe practices during the COVID 19 pandemic. Local conditions will almost certainly dictate modifications. Hopefully there is sufficient information to empower caregivers to make good decisions to keep patients, co-workers and themselves safe.

Breathing Circuit Filter Virus-Safety Product and Manufacturer Information

<https://www.apsf.org/faq-on-anesthesia-machine-use-protection-and-decontamination-during-the-covid-19-pandemic/>

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





Aespire 7900

User's Reference Manual—Part 2

Setup, Cleaning and Sterilization, Maintenance and Troubleshooting


Breathing system cleanable parts

- 
1. Bellows assembly
 2. APL Valve Ramp
 3. APL Diaphragm
 4. Breathing circuit module (O₂ sensor not autoclavable)
 5. Absorber canister (reusable only)
 6. Flow sensor cover
 7. Flow Sensor Module (plastic flow sensors not autoclavable)
 8. Exhalation Valve assembly

Figure 2-1 • Autoclavable assemblies

Parts marked 134°C are autoclavable or washable by hand or machine (mild detergent pH <10.5). Rinse and dry completely. All parts except the O₂ sensor and disposable flow sensors can be washed.

If the flow sensors are plastic, refer to the "How to clean and disinfect the flow sensors" procedure. If metal, they can be autoclaved at 134°C.

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